

Immunocapture Real-Time PCR to Detect Mycotoxigenic Mold Spores in Grains and Foods


M. A. Cousin
Department of Food Science
C. P. Woloshuk
Department of Botany and Plant Pathology
Purdue University

Mycotoxigenic Molds

- *Aspergillus* species: aflatoxins (B₁ and M₁), ochratoxin A
- *Fusarium* species: deoxynivalenol (DON), fumonisins (B₁ and B₂), zearalenone
- *Penicillium* species: ochratoxin A
- Action levels for aflatoxins
- Advisory levels for DON, fumonisins, zearalenone

Detect Molds or Mycotoxins?

Mycotoxins can be produced during some food manufacture.



Previous Research

- Develop immunocapture real-time PCR method for *Fusarium* species
- Limitation on PCR due to inability to release DNA from spores

Objective of Research

- To study methods to break mold spores to release DNA
- To incorporate the best method into immunocapture real-time PCR to detect *molds* in foods and grains
- To develop a PCR library to mycotoxigenic molds

Methods to Study for Spore Breakage

- Enzymatic - lyticase digestion
- Physical methods - heat or cold
- Mechanical methods – ball mill, FastPrep™
- Partial spore germination
- Combination of methods

Breakage of Spores

- Preliminary research
 1. Lyticase can break some *Fusarium* spores
 2. Some *Fusarium* spores germinate in 2-3 h
 3. FastPrep™ variable for *Fusarium* species
- Physical methods – just beginning

PCR Library For Detection of Mycotoxigenic Molds

- Current Activities
 - Multiplex Genus-specific Assay: Developing TaqMan primers and probes to detect the presence of *Penicillium*, *Fusarium*, and *Aspergillus* species
 - Multiplex Toxin-specific Assay: Assembling TaqMan primers and probes to detect pathway genes for DON, fumonisin, ochratoxin, and aflatoxin

Future Research

- Incorporate procedures to break spores into real-time PCR
- Simplify immunocapture before real-time PCR
- Test in industry to determine cost effectiveness